



Pulverizing and Blending

SPEX SamplePrep 8530 Enclosed Shatterbox®

8530-230 Enclosed Shatterbox

For 230 V/50 Hz operation,
CE Approved.

- Features a cam-activated clamp that requires only moderate hand pressure to secure the grinding dish in the Shatterbox and is adjustable to secure containers of different heights.
- LCD display with environmentally safe, push-button membrane switch.
- Soundproof enclosure reduces noise levels.
- Equipped with a safety interlock system and lockable casters.
- Cooling fan to prevent the motor from overheating.
- Wide choice of grinding containers, sold separately.
- New model 8530-115 for 115 V/60 Hz operation coming soon.



8530 Enclosed Shatterbox

Specifications:

| | |
|-------------|---|
| Dimensions: | 19 in. (48 cm) x 19 in. (48 cm) x 40 in. (102 cm) |
| Weight: | 225 lbs. net (102 kg); 246 lbs. gross (112 kg) |
| Voltage: | 230 V/50 Hz |
| Motor: | ½ horsepower @ 835 RPM |