















LISTA DE PRECIOS KERN 2010

Modelo	Descripción	Rango peso (g)	Resolución (g)	Dimensiones del plato (mm)	Precio 2010	
Precision Balances						
	440-21A	Precision balance	60	0,00100	81	290,00 €
	440-21N	Precision balance	40	0,00100	81	290,00 €
	440-33N	Precision balance	200	0,01000	105	240,00 €
	440-35A	Precision balance	600	0,01000	105	300,00 €
	440-35N	Precision balance	400	0,01000	105	270,00 €
	440-43N	Precision balance	400	0,10000	130 x 130	150,00 €
	440-45N	Precision balance	1.000	0,10000	130 x 130	200,00 €
	440-47N	Precision balance	2.000	0,10000	130 x 130	220,00 €
	440-49A	Precision balance	6.000	0,10000	150 x 170	300,00 €
	440-49N	Precision balance	4.000	0,10000	150 x 170	270,00 €
	440-51N	Precision balance	4.000	1,00000	150 x 170	130,00 €
	440-53N	Precision balance	6.000	1,00000	150 x 170	150,00 €
	440-55N	Precision balance	6.000	0,20000	150 x 170	240,00 €
	440-210-002	Protective working cover for für 440-21N	-	-	-	20,00 €
	440-450-002	Protective working cover for 440-43N,440-45N,440-47N,KB 1200-2, KB 1600-2	-	-	-	20,00 €
	440-530-002	Protective Working cover for 440-49N,440-49A,440-51N,440-53N,440-55N,KB 10000-1	-	-	-	20,00 €
	440-902	Adaptor EU for 440, 442, DE	-	-	-	25,00 €
	440-A01	Hook for underfloor weighing	-	-	-	20,00 €
	470-902-002	Mains adaptor UK for series EW/EG, 440, DE, CB	-	-	-	25,00 €
	470-902-003	Mains adaptor US 110V/60Hz for series 440, CB, DE, ECB, EMB, EOB	-	-	-	25,00 €
	572-30	Precision balance	161	0,00100	105	380,00 €
	572-31	Precision balance	241	0,00100	105	390,00 €
	572-32	Precision balance	421	0,00100	105	420,00 €
	572-33	Precision balance	1.610	0,01000	150	360,00 €
	572-35	Precision balance	2.410	0,01000	150	370,00 €
	572-37	Precision balance	3.010	0,01000	150	380,00 €
	572-39	Precision balance	4.210	0,01000	150	430,00 €
	572-43	Precision balance	10.100	0,10000	160 x 200	350,00 €
	572-45	Precision balance	12.100	0,05000	160 x 200	380,00 €
	572-49	Precision balance	16.100	0,10000	160 x 200	360,00 €
	572-55	Precision balance	20.100	0,05000	160 x 200	430,00 €
	572-57	Precision balance	24.100	0,10000	160 x 200	380,00 €
	572-926	Interface cable, balance - PC FCB, PCB, NDE, FKB-A, 440, 572, 573, CB, DE, DS, KB, ITB - ITT, FTB, FTC	-	-	-	46,00 €
	572-A02	Protective working cover over key pad and casing	-	-	-	25,00 €
572-A05	Glass draught guard with sliding doors for KERN 572-31	-	-	-	160,00 €	
	573-34NM	Precision balance, with type approval	650	0,01000	160 x 200	450,00 €
	573-36NM	Precision balance, with type approval	2.000	0,01000	160 x 200	480,00 €
	573-46NM	Precision balance, with type approval	6.500	0,10000	160 x 200	450,00 €
	573-56NM	Precision balance, with type approval	20.100	0,10000	160 x 200	480,00 €
	EG 2200-2NM	Precision balance, with type approval	2.200	0,01000	180 x 160	740,00 €
EG 220-3NM	Precision balance, with type approval	220	0,00100	118	720,00 €	
EG 4200-2NM	Precision balance, with type approval	4.200	0,01000	180 x 160	830,00 €	
EG 420-3NM	Precision balance, with type approval	420	0,00100	118	810,00 €	
EG 620-3NM	Precision balance, with type approval	620	0,00100	118	930,00 €	
EG-A01	Full glass draught guard with sliding doors	-	-	-	160,00 €	
EG-A02	Protective working cover, standard for series EW/EG	-	-	-	25,00 €	
EG-A03	Glass draft shield only for EG ...-3NM	-	-	-	160,00 €	
EG-A04	Rechargeable battery pack, internal series EG-N and EW-N < Max. 820g	-	-	-	100,00 €	
EG-A05	Protective working cover for models EG-NM	-	-	-	25,00 €	
EG-A06	Rechargeable battery pack, internal for series EG-N, EW-N, Max > 2200g	-	-	-	100,00 €	
EG-A07	Hook for underfloor weighing	-	-	-	36,00 €	
EG-A08	Hook for underfloor weighing	-	-	-	36,00 €	
EG-A09	Protective working cover, standard for series EW/EG	-	-	-	25,00 €	
	EW 1500-2M	Precision balance, with type approval	1.500	0,01000	140	570,00 €
	EW 150-3M	Precision balance, with type approval	150	0,00100	110	560,00 €
	EW 3000-2M	Precision balance, with type approval	3.000	0,01000	140	670,00 €
	EW 6000-1M	Precision balance, with type approval	6.000	0,10000	172 x 142	570,00 €
	EW 600-2M	Precision balance, with type approval	600	0,01000	140	520,00 €
	EW-A01	Data output RS 232 for series EW	-	-	-	120,00 €
	EW-A02	Rechargeable battery pack, internal for series EW	-	-	-	100,00 €
	EW-A02-001	Accu-kit for series EW	-	-	-	34,36 €
	EW-A05	Carat pan approx. 300 g, diameter 290mm	-	-	-	50,00 €
	474-926	Interface cable balance-PC for series 474, EW, EG, EC	-	-	-	46,00 €
	EW 12000-1NM	Precision balance, with type approval	12.000	0,10000	180 x 160	670,00 €
	EW 220-3NM	Precision balance, with type approval	600	0,01000	118	610,00 €
	EW 2200-2NM	Precision balance, with type approval	2.200	0,01000	180 x 160	640,00 €
	EW 4200-2NM	Precision balance, with type approval	4.200	0,01000	180 x 160	720,00 €
	EW 420-3NM	Precision balance, with type approval	420	0,00100	118	680,00 €
	EW 6200-2NM	Precision balance, with type approval	6.200	0,01000	180 x 160	870,00 €
	EW 620-3NM	Precision balance, with type approval	620	0,00100	118	800,00 €
	EW 820-2NM	Precision balance, with type approval	820	0,01000	170 x 142	560,00 €








LISTA DE PRECIOS KERN 2010

	Modelo	Descripción	Rango peso (g)	Resolución (g)	Dimensiones del plato (mm)	Precio 2010
	EWB 1200-1M	Precision balance, with type approval	1.200	0,10000	170 x 140	440,00 €
	EWB 220-2M	Precision balance, with type approval	220	0,01000	140	440,00 €
	EWB 620-2M	Precision balance, with type approval	620	0,01000	140	470,00 €
	EWB-A01	Rechargeable battery pack, internal for series EWB	-	-	-	100,00 €
	KB 10000-1N	Precision balance	10.100	0,10000	150 x 170	230,00 €
	KB 10K0.05N	Precision balance	10.100	0,05000	150 x 170	260,00 €
	KB 10000-1NM	Precision balance, with type approval	10.000	0,10000	150 x 170	380,00 €
	KB 1200-2N	Precision balance	1.210	0,01000	130 x 130	240,00 €
	KB 120-3N	Precision balance	121	0,00100	81	240,00 €
	KB 2000-2N	Precision balance	2.010	0,01000	130 x 130	260,00 €
	KB 2000-2NM	Precision balance, with type approval	2.010	0,01000	130 x 130	380,00 €
	KB 200-3N	Precision balance	201	0,00100	81	260,00 €
	KB 2400-2N	Precision balance	2.410	0,01000	130 x 130	270,00 €
	KB 240-3N	Precision balance	241	0,00100	81	270,00 €
	KB 3600-2N	Precision balance	3.600	0,01000	130 x 130	290,00 €
	KB 360-3N	Precision balance	361	0,00100	81	290,00 €
	KB 6500-1NM	Precision balance, with type approval	6.500	0,10000	150 x 170	340,00 €
	KB 650-2NM	Precision balance, with type approval	650	0,01000	130 x 130	340,00 €
	PBJ 4200-2M	Precision balance, with type approval	4.200	0,01000	180 x 170	1.110,00 €
	PBJ 420-3M	Precision balance, with type approval	420	0,00100	105 x 105	1.050,00 €
	PBJ 6200-2M	Precision balance	6.200	0,01000	180 x 170	1.280,00 €
	PBJ 620-3M	Precision balance, with type approval	620	0,00100	105 x 105	1.190,00 €
	PBJ 8200-1M	Precision balance, with type approval	8.200	0,01000	180 x 170	1.050,00 €
	PBS 4200-2M	Precision balance	4.200	0,01000	180 x 170	950,00 €
	PBS 420-3M	Precision balance	420	0,00100	105 x 105	890,00 €
	PBS 6200-2M	Precision balance	6.200	0,01000	180 x 170	1.080,00 €
	PBS 620-3M	Precision balance	620	0,00100	105 x 105	1.020,00 €
PBS 8200-1M	Precision balance	8.200	0,10000	180 x 170	920,00 €	
	PCB 100-3	Precision balance	100	0,00100	81	210,00 €
	PCB 200-2	Precision balance	200	0,01000	105	170,00 €
	PCB 250-3	Precision balance	250	0,00100	81	240,00 €
	PCB 400-1	Precision balance	400	0,10000	130 x 130	130,00 €
	PCB 400-2	Precision balance	400	0,01000	105	150,00 €
	PCB 350-3	Precision balance	350	0,00100	81	270,00 €
	PCB 1000-1	Precision balance	1.000	0,10000	130 x 130	135,00 €
	PCB 1000-2	Precision balance	1.000	0,01000	130 x 130	195,00 €
	PCB 2000-1	Precision balance	2.000	0,10000	130 x 130	160,00 €
	PCB 2500-2	Precision balance	2.500	0,01000	130 x 130	240,00 €
	PCB 3500-2	Precision balance	3.500	0,01000	130 x 130	260,00 €
	PCB 4000-1	Precision balance	4.000	0,10000	150 x 170	190,00 €
	PCB 6000-0	Precision balance	6.000	1,00000	150 x 170	100,00 €
	PCB 6000-1	Precision balance	6.000	0,10000	150 x 170	180,00 €
	PCB 10000-1	Precision balance	10.000	0,10000	150 x 170	190,00 €
	PCB 10K0.05	Precision balance	10.000	0,05000	150 x 170	220,00 €
	PCB-A01	Internal rechargeable battery pack EURO for series KERN PCB, FKB-A, ECB, ECE	-	-	-	45,00 €
	PCB-A02	Protective Working cover for PCB with platform Ø 80 mm	-	-	-	20,00 €
	PCB-A03	Protective Working cover for PCB with platform Ø 105 mm	-	-	-	20,00 €
	PCB-A04	Protective Working cover for PCB with platform 130 x 130 mm	-	-	-	20,00 €
	PCB-A05	Protective Working cover for PCB with platform 150 x 170 mm	-	-	-	20,00 €
	PCB-A06	Rechargeable battery pack UK, internal	-	-	-	50,00 €
	PCB-A07B	Rechargeable battery pack, internal for series PCB-B	-	-	-	50,00 €






LISTA DE PRECIOS KERN 2010

	Modelo	Descripción	Rango peso (g)	Resolución (g)	Dimensiones del plato (mm)	Precio 2010	
	PEJ 2200-2M	Precision balance, with type approval	2.200	0,01000	200 x 200	1.100,00 €	
	PEJ 220-3M	Precision balance, with type approval	220	0,00100	140 x 120	930,00 €	
	PEJ 4200-2M	Precision balance, with type approval	4.200	0,01000	200 x 200	1.220,00 €	
	PEJ 420-3M	Precision balance, with type approval	420	0,00100	140 x 120	1.060,00 €	
	PEJ 620-3M	Precision balance, with type approval	620	0,00100	140 x 120	1.300,00 €	
	PES 15000-1M	Precision balance, with type approval	15.000	0,10000	200 x 200	1.250,00 €	
	PES 2200-2M	Precision balance, with type approval	2.200	0,01000	200 x 200	890,00 €	
	PES 220-3M	Precision balance, with type approval	220	0,00100	140 x 120	740,00 €	
	PES 31000-1M	Precision balance, with type approval	31.000	0,10000	250 x 220	1.990,00 €	
	PES 4200-2M	Precision balance, with type approval	4.200	0,01000	200 x 200	1.050,00 €	
	PES 420-3M	Precision balance, with type approval	420	0,00100	140 x 120	870,00 €	
	PES 6200-2M	Precision balance, with type approval	6.200	0,01000	200 x 200	1.250,00 €	
	PES 620-3M	Precision balance, with type approval	620	0,00100	140 x 120	1.090,00 €	
	PES 8200-1M	Precision balance, with type approval	8.200	0,10000	200 x 200	1.150,00 €	
	PES-A01	Internal rechargeable battery pack EURO for series KERN PES/PEJ	-	-	-	-	560,00 €
	PFB 120-3	Precision balance	120	0,00100	80	210,00 €	
	PFB 200-3	Precision balance	200	0,00100	80	235,00 €	
	PFB 1200-2	Precision balance	1.200	0,01000	120	180,00 €	
	PFB 2000-2	Precision balance	2.000	0,01000	120	235,00 €	
	PFB 6K0.05	Precision balance	6.000	0,05000	155 x 145	235,00 €	
	PFB 6000-1	Precision balance	6.000	0,10000	155 x 145	170,00 €	
	PLE 310-3N	Precision balance	310	0,00100	80	320,00 €	
	PLE 3100-2N	Precision balance	3.100	0,01000	135	320,00 €	
	PLS 310-3F	Precision balance	310	0,00100	110	410,00 €	
	PLS 510-3A	Precision balance	510	0,00100	110	690,00 €	
	PLS 720-3A	Precision balance	720	0,00100	110	730,00 €	
	PLS 1200-3A	Precision balance	1.200	0,00100	110	790,00 €	
	PLS 1200-3DA	Precision balance	720; 1.200	0,00100;	110	700,00 €	
	PLS 3100-2F	Precision balance	3.100	0,01000	160	410,00 €	
	PLS 4200-2A	Precision balance	4.200	0,01000	160	690,00 €	
	PLS 6200-2A	Precision balance	6.200	0,01000	160	730,00 €	
	PLS 6200-2DA	Precision balance	4.200; 6.200	0,01000;	160	700,00 €	
	PLJ 360-3NM	Precision balance, with type approval	360	0,00100	128 x 128	750,00 €	
	PLJ 600-3NM	Precision balance, with type approval	600	0,00100	128 x 128	810,00 €	
	PLJ 600-2GM	Precision balance, with type approval	600	0,01000	128 x 128	590,00 €	
	PLJ 750-3N	Precision balance, with type approval	750	0,00100	128 x 128	990,00 €	
	PLJ 3500-2NM	Precision balance, with type approval	3.500	0,01000	195 x 195	750,00 €	
	PLJ 4500-2NM	Precision balance, with type approval	4.500	0,01000	195 x 195	810,00 €	
	PLJ 6000-1GM	Precision balance, with type approval	6.000	0,10000	195 x 195	590,00 €	
	PLJ-A01	Protective working cover	-	-	-	-	25,00 €
	PLJ-A02	Hook for suspended weighing for series KERN PL-A/PL-F	-	-	-	-	20,00 €
	PLJ-A03	Mains adaptor EURO for KERN series PLS/PLJ	-	-	-	-	35,00 €
	PLJ-A04	Mains adaptor UK for KERN series PLS/PLJ	-	-	-	-	35,00 €
PLJ-A05	Mains adaptor US for KERN series PLS/PLJ	-	-	-	-	35,00 €	
	PRJ 10200-1NM	Precision balance, with type approval	10.200	0,10000	200 x 200	1.780,00 €	
	PRJ 10200-1IP65NM	Precision balance, with type approval	10.200	0,10000	200 x 200	1.930,00 €	
	PRJ 1200-3N	Precision balance	1.200	0,00100	135 x 135	2.380,00 €	
	PRJ 320-3NM	Precision balance, with type approval	320	0,00100	135 x 135	1.850,00 €	
	PRJ 4200-2NM	Precision balance, with type approval	4.200	0,01000	170 x 170	1.870,00 €	
	PRJ 4200-2IP65NM	Precision balance, with type approval	4.200	0,01000	170 x 170	2.020,00 €	
	PRJ 6200-2NM	Precision balance, with type approval	6.200	0,01000	170 x 170	2.100,00 €	
	PRJ 620-3NM	Precision balance, with type approval	620	0,00100	135 x 135	2.040,00 €	
	PRJ 8200-1IP65NM	Precision balance	8.200	0,10000	200 x 200	1.850,00 €	
	PRJ 8200-1NM	Precision balance, with type approval	8.200	0,10000	200 x 200	1.700,00 €	
	PRS 12200-1N	Precision balance	12.200	0,10000	200 x 200	1.440,00 €	
	PRS 12200-1IP65N	Precision balance	12.200	0,10000	200 x 200	1.590,00 €	
	PRS 320-3N	Precision balance	320	0,00100	135 x 135	1.550,00 €	
	PRS 4200-2N	Precision balance	4.200	0,01000	170 x 170	1.550,00 €	
	PRS 4200-2IP65N	Precision balance	4.200	0,01000	170 x 170	1.700,00 €	
	PRS 6200-2N	Precision balance	6.200	0,01000	170 x 170	1.790,00 €	
	PRS 620-3N	Precision balance	620	0,00100	135 x 135	1.730,00 €	
	PRS 8200-1N	Precision balance	8.200	0,10000	200 x 200	1.400,00 €	
	PRS 8200-1IP65N	Precision balance	8.200	0,10000	200 x 200	1.550,00 €	
	PR-A01	Draftshield 180 mm for series PRS/PRJ	-	-	-	-	215,00 €
	PR-A04	Hook for suspended weighing	-	-	-	-	45,00 €
	PR-A05	Suppression holder	-	-	-	-	95,00 €
	PR-A06	Adapter RJ45/DB25 Male (Printer)	-	-	-	-	75,00 €
	PR-A07	Adapter RJ45/DB9 Female (PC)	-	-	-	-	75,00 €
	PR-A08	-	-	-	-	-	75,00 €
	PR-A09	Remote display for the table with cable 1,5 m	-	-	-	-	395,00 €
	PR-A10	Second remote display inclusive stand	-	-	-	-	430,00 €
	PR-A11	Data interface, Second RS 232 C	-	-	-	-	230,00 €
	PR-A12	Data interface 20 mA current loop	-	-	-	-	230,00 €
	PR-A13	Data interface, Analogue, -10 V ... +10 V	-	-	-	-	370,00 €
	PR-A14	Smartbox for parts counting	-	-	-	-	530,00 €
	PR-A15	Smartbox for prepackaging	-	-	-	-	530,00 €
PR-A16	CCD-barcode for Smart-Box (KERN PRS/PRJ)	-	-	-	-	1.210,00 €	
PR-A18	Indicator light for smartbox	-	-	-	-	930,00 €	
PR-A19	Alarm horn for indicator light	-	-	-	-	185,00 €	
PR-A20	-	-	-	-	-	335,00 €	
PR-A21	Multiplexer for 7 balances (RS 232)	-	-	-	-	360,00 €	
PR-A22	Accu set, complete, rechargeable	-	-	-	-	295,00 €	
PR-A23	Interface cable for PC	-	-	-	-	130,00 €	
PR-A24	Interface cable for printer	-	-	-	-	130,00 €	



LISTA DE PRECIOS KERN 2010

	Modelo	Descripción	Rango peso (g)	Resolución (g)	Dimensiones del plato (mm)	Precio 2010
	PR-A25	Protective working cover for series PRS/PRJ	-	-	-	12,00 €
	PR-A31	Relay box	-	-	-	330,00 €
	EMB 100-3	School balances	100	0,00100	82	150,00 €
	EMB 200-2	School balances	200	0,01000	105	100,00 €
	EMB 200-3	School balances	200	0,00100	82	170,00 €
	EMB 220-1	School balances	220	0,10000	150	40,00 €
	EMB 500-1	School balances	500	0,10000	150	44,00 €
	EMB 600-2	School balances	600	0,01000	105	125,00 €
	EMB 1000-2	School balances	1.000	0,01000	150	150,00 €
	EMB 1200-1	School balances	1.200	0,10000	150	80,00 €
	EMB 2000-2	School balances	2.000	0,01000	105	170,00 €
	EMB 2200-0	School balances	2.200	1,00000	150	40,00 €
	EMB 5.2K1	School balances	5.200	1,00000	150	50,00 €
	EMB 5.2K5	School balances	5.200	5,00000	150	40,00 €
	EMB 6000-1	School balances	6.000	0,10000	150	125,00 €
	EMB-A02	Stainless steel weighing pan 105 mm	-	-	-	25,00 €
EMB-A03	Stainless steel weighing pan 150 mm	-	-	-	25,00 €	
Analytical Balances						
	ABJ 120-4M	Analytical balance, with type approval	120	0,00010	80	880,00 €
	ABJ 220-4M	Analytical balance, with type approval	220	0,00010	80	960,00 €
	ABJ 320-4	Analytical balance, with type approval	320	0,00010	80	1.190,00 €
	ABJ 80-4M	Analytical balance, with type approval	83	0,00010	80	830,00 €
	ABS 120-4	Analytical balance	120	0,00010	80	780,00 €
	ABS 220-4	Analytical balance	220	0,00010	80	860,00 €
	ABS 80-4	Analytical balance	83	0,00010	80	730,00 €
ABS-A02	Set for the determination of density of liquids and solids (models with d =0,001g)	-	-	-	650,00 €	
ABS-A05	Interface cable for PC for models ABS	-	-	-	80,00 €	
	ABT 120-4M	Analytical balance, with type approval	120	0,00010	80	1.390,00 €
	ABT 120-5DM	Analytical balance, with type approval	42; 120	0,00001;	80	1.750,00 €
	ABT 100-5M	Analytical balance, with type approval	101	0,00001	80	1.690,00 €
	ABT 220-4M	Analytical balance, with type approval	220	0,00010	80	1.510,00 €
	ABT 220-5DM	Analytical balance, with type approval	82; 220	0,00001;	80	1.950,00 €
	ABT 320-4M	Analytical balance, with type approval	320	0,00010	80	1.690,00 €
	ABT-A01	Set for the determination of density of liquids and solids (models with d =0,001g)	-	-	-	640,00 €
ABT-A02	Protective working cover for models ABT	-	-	-	25,00 €	
ABT-P01	Package Solution Pipettes Calibration EASYCAL	-	-	-	2.350,00 €	
	AEJ 120-4M	Analytical balance, with type approval	120	0,00010	80	1.240,00 €
	AEJ 220-4M	Analytical balance, with type approval	220	0,00010	80	1.300,00 €
	AES 120-4	Analytical balance	120	0,00010	80	1.170,00 €
	AES 220-4	Analytical balance	220	0,00010	80	1.220,00 €
	AES-A01	Density determination kit for series KERN AES/AEJ	-	-	-	460,00 €
Moisture Analyzers						
	MRS 120-3	Moisture analyzer	120	0,00100	-	2.620,00 €
	MRS-A01	Temperature calibration set consisting of measuring sensor and analysis box	-	-	-	630,00 €
	MLS 65-3A	Moisture analyzer	65	0,00100	-	1.190,00 €
	MLS 50-3HA160N	Moisture analyzer	50	0,00100	-	1.550,00 €
	MLB 50-3N	Moisture analyzer	50	0,00100	-	1.220,00 €
	MLB 60-3A	Moisture analyzer	60	0,00100	-	1.220,00 €
	MLB-A01	Plates for samples, of Aluminium, 80 pieces, f. moist.analyzer MLB	-	-	-	45,00 €
	MLB-A05	Interface cable, balance - PC, for series ALS/ALJ, PLS/PLJ, MLB, MLS, ALT	-	-	-	46,00 €
MLB-A11	Temperature calibration set, consists of measurement sensor and display device	-	-	-	250,00 €	

LISTA DE PRECIOS KERN 2010

	Modelo	Descripción	Rango peso (g)	Resolución (g)	Dimensiones del plato (mm)	Precio 2010
Bench Scales						
	FCB 15K1BDM	Bench scale	7.500;	1,00000;	253 x 228	520,00 €
	FCB 8K0.1BM	Bench scale	8.000	0,10000	253 x 228	510,00 €
	FCB 15K0.1BM	Bench scale	15.000	0,10000	253 x 228	540,00 €
	FCB 12K0.05B	Bench scale	12.100	0,05000	252 x 228	370,00 €
	FCB 12K0.1B	Bench scale	12.100	0,10000	252 x 228	340,00 €
	FCB 12K1	Bench scale	12.000	1,00000	252 x 228	150,00 €
	FCB 15K5	Bench scale	15.000	5,00000	252 x 228	130,00 €
	FCB 16K0.2	Bench scale	16.000	0,20000	252 x 228	195,00 €
	FCB 24K0.1B	Bench scale	24.100	0,10000	252 x 228	370,00 €
	FCB 24K0.2B	Bench scale	24.100	0,20000	252 x 228	340,00 €
	FCB 24K10	Bench scale	24.000	10,00000	252 x 228	130,00 €
	FCB 24K2	Bench scale	24.000	2,00000	252 x 228	150,00 €
	FCB 30K0.5	Bench scale	30.000	0,50000	252 x 228	195,00 €
	FCB 30K1	Bench scale	30.000	1,00000	252 x 228	185,00 €
	FCB 3K0.1	Bench scale	3.000	0,10000	252 x 228	185,00 €
	FCB 6K0.02B	Bench scale	6.100	0,02000	252 x 228	390,00 €
	FCB 6K0.5	Bench scale	6.000	0,50000	252 x 228	160,00 €
	FCB 8K0.1	Bench scale	8.000	0,10000	252 x 228	210,00 €
	FCB-A01	Rechargeable battery pack, internal for series 440, CB, DE, FCB	-	-	-	35,00 €
	FCB-A02	Protective working cover for series FCB, CB, DE	-	-	-	20,00 €
FCB-A03	Rechargeable battery pack UK, internal for series 440, CB, DE, FCB	-	-	-	35,00 €	
FCB-A03B	Rechargeable battery pack, internal for series FCB-B	-	-	-	40,00 €	
	FCE 15K5	Bench scale	15.000	5,00000	252 x 228	105,00 €
	FCE 30K10	Bench scale	30.000	10,00000	252 x 228	105,00 €
	FCE 30K10L	Bench scale	30.000	10,00000	340 x 240	150,00 €
	FCE 3K1	Bench scale	3.000	1,00000	252 x 228	105,00 €
	FCE 60K20	Bench scale	60.000	20,00000	340 x 240	150,00 €
	FCE 6K2	Bench scale	6.000	2,00000	252 x 228	105,00 €
	FFN 3K0.5IP	Bench scale	3.000	0,50000	230 x 190	260,00 €
	FFN 6K1IP	Bench scale	6.000	1,00000	231 x 190	260,00 €
	FFN 15K2IP	Bench scale	15.000	2,00000	232 x 190	260,00 €
	FFN 25K5IP	Bench scale	25.000	5,00000	233 x 190	260,00 €
	FKB 15K0.1M	Bench scale	15.000	0,10000	340 x 240	580,00 €
	FKB 15K0.5A	Bench scale	15.000	0,50000	340 x 240	240,00 €
	FKB 15K1A	Bench scale	15.000	1,00000	340 x 240	210,00 €
	FKB 16K0.05	Bench scale	16.100	0,05000	340 x 240	390,00 €
	FKB 16K0.1	Bench scale	16.100	0,10000	340 x 240	370,00 €
	FKB 16K0.2A	Bench scale	16.000	0,20000	340 x 240	250,00 €
	FKB 30K1A	Bench scale	30.000	1,00000	340 x 240	230,00 €
	FKB 30K2A	Bench scale	30.000	2,00000	340 x 240	210,00 €
	FKB 30K5DA	Bench scale	15.000;	5,00000;	340 x 240	190,00 €
			30.000	10,00000		
	FKB 36K0.1	Bench scale	36.100	0,10000	340 x 240	420,00 €
	FKB 36K0.2	Bench scale	36.100	0,20000	340 x 240	390,00 €
	FKB 36K0.5A	Bench scale	36.000	0,50000	340 x 240	250,00 €
	FKB 50K1TM	Bench scale	10.000	1,00000	340 x 240	580,00 €
	FKB 60K10DA	Bench scale	30.000;	10,00000;	340 x 240	195,00 €
			60.000	20,00000		
	FKB 60K2A	Bench scale	60.000	2,00000	340 x 240	255,00 €
	FKB 60K5A	Bench scale	60.000	5,00000	340 x 240	210,00 €
	FKB 65K0.2	Bench scale	65.100	0,20000	340 x 240	420,00 €
	FKB 65K0.5	Bench scale	65.100	0,50000	340 x 240	390,00 €
	FKB 65K1A	Bench scale	65.000	1,00000	340 x 240	260,00 €
	FKB 65K1M	Bench scale	65.000	1,00000	340 x 240	580,00 €
FKB 6K0.02	Bench scale	6.100	0,02000	340 x 240	390,00 €	
FKB 8K0.05	Bench scale	8.100	0,05000	340 x 240	370,00 €	
FKB 8K0.1M	Bench scale	8.000	0,10000	340 x 240	530,00 €	
FKB 8K0.1A	Bench scale	8.000	0,10000	340 x 240	250,00 €	
FKB-A02	Protective working cover for series FKB; CKE; NKE; OKE	-	-	-	20,00 €	
	FKT 6K0.02	Bench scale	6.100	0,02000	340 x 240	550,00 €
	FKT 6K0.02L	Bench scale	6.100	0,02000	340 x 240	590,00 €
	FKT 12K0.05	Bench scale	12.100	0,05000	253 x 228	550,00 €
	FKT 16K0.05L	Bench scale	16.100	0,05000	340 x 240	590,00 €
	FKT 24K0.1	Bench scale	24.100	0,10000	253 x 228	550,00 €
	FKT 36K0.1L	Bench scale	36.100	0,10000	340 x 240	610,00 €
FKT 65K0.2L	Bench scale	65.100	0,20000	340 x 240	610,00 €	

LISTA DE PRECIOS KERN 2010

	Modelo	Descripción	Rango peso (g)	Resolución (g)	Dimensiones del plato (mm)	Precio 2010
	FTB 15K0.5	Bench scale	15.000	0,50000	240 x 200	670,00 €
	FTB 15K0.5IP	Bench scale	15.000	0,50000	240 x 200	890,00 €
	FTB 15K0.5L	Bench scale	15.000	0,50000	350 x 240	770,00 €
	FTB 15K2DLM	Bench scale	6.000;	2,00000;	350 x 240	770,00 €
	FTB 15K2DM	Bench scale	6.000;	2,00000;	240 x 200	670,00 €
	FTB 35K1	Bench scale	35.000	1,00000	240 x 200	770,00 €
	FTB 35K5DM	Bench scale	15.000;	5,00000;	350 x 240	770,00 €
	FTB 3K0.1	Bench scale	3.000	0,10000	240 x 200	670,00 €
	FTB 3K0.5DM	Bench scale, with type approval	1.500; 3.000	0,50000;	240 x 200	670,00 €
	FTB 60K10DM	Bench scale, with type approval	30.000;	10,00000;	350 x 240	770,00 €
	FTB 60K2	Bench scale	60.000	2,00000	350 x 240	770,00 €
	FTB 6K0.2	Bench scale	6.000	0,20000	240 x 200	670,00 €
	FTB 6K1DM	Bench scale, with type approval	3.000; 6.000	1,00000;	240 x 200	670,00 €
	FTB 15K0.5LIP	Bench scale	15.000	0,50000	350 x 240	1.000,00 €
	FTB 15K2DLMIP	Bench scale, with type approval	15.000	2,00000	350 x 240	1.000,00 €
	FTB 15K2DMIP	Bench scale, with type approval	6.000;	2,00000;	200 x 240	900,00 €
	FTB 35K1IP	Bench scale	35.000	1,00000	350 x 240	1.000,00 €
	FTB 35K5DMIP	Bench scale, with type approval	15.000;	5,00000;	350 x 240	1.000,00 €
	FTB 3K0.1IP	Bench scale	3.000	1,00000	200 x 240	900,00 €
	FTB 3K0.5DMIP	Bench scale, with type approval	1.500; 3.000	0,50000;	200 x 240	900,00 €
	FTB 60K10DMIP	Bench scale, with type approval	30.000;	10,00000;	350 x 240	1.000,00 €
	FTB 60K2IP	Bench scale	60.000	2,00000	350 x 240	1.000,00 €
	FTB 6K0.2IP	Bench scale	6.000	2,00000	200 x 240	900,00 €
	FTB 6K1DMIP	Bench scale, with type approval	6.000	1,00000	200 x 240	900,00 €
	FTB-A01-001	Rechargeable battery pack, internal factory option, for weighing plates 240x200mm	-	-	-	200,00 €
	FTB-A01-002	Rechargeable battery pack, internal for weighing plates 240x200mm	-	-	-	250,00 €
	FTB-A02-001	Protective working cover (order with balance) for weighing plates 240x200mm	-	-	-	38,00 €
	FTB-A02-002	Protective working cover (for reorders) for weighing plates 240x200mm	-	-	-	40,00 €
	FTB-A03-001	Protective working cover (order with balance) for weighing plates 350x240mm	-	-	-	40,00 €
	FTB-A03-002	Protective working cover (for reorders) for weighing plates 350x240mm	-	-	-	45,00 €
	FTB-A04-001		-	-	-	200,00 €
	FTB-A04-002		-	-	-	250,00 €
FTB-A05-001	Rechargeable battery pack, internal factory option, for weighing plates 240x200mm	-	-	-	320,00 €	
FTB-A05-002	Rechargeable battery pack, internal for weighing plates 240x200mm	-	-	-	320,00 €	
FTB-A06-001	Rechargeable battery pack, internal factory option, for "large" housing/casing	-	-	-	320,00 €	
FTB-A06-002	Rechargeable battery pack, internal factory option, for "large" housing/casing	-	-	-	320,00 €	
FTB-A09	Ethernet interface for series KERN FTB	-	-	-	410,00 €	
	GAB 3K1DM	Bench scale, with type approval	1.500; 3.000	1,00000;	230 x 180	160,00 €
	GAB 6K2DM	Bench scale, with type approval	3.000; 6.000	2,00000;	230 x 180	160,00 €
	GAB 6K1DNM	Bench scale, with type approval	3.000; 6.000	1,00000;	294 x 225	210,00 €
	GAB 6K0.05N	Bench scale, with type approval	6.000	0,05000	294 x 225	250,00 €
	GAB 15K2DNM	Bench scale, with type approval	6.000;	2,00000;	294 x 225	210,00 €
	GAB 15K5DM	Bench scale, with type approval	15.000	5,00000	280 x 260	170,00 €
	GAB 12K0.1N	Bench scale, with type approval	12.000	0,10000	294 x 225	250,00 €
	GAB 30K10DM	Bench scale, with type approval	15.000;	10,00000;	280 x 260	170,00 €
	GAB 30K5DNM	Bench scale, with type approval	30.000	20,00000		
	GAB 30K5DNM	Bench scale, with type approval	15.000;	5,00000;	294 x 225	210,00 €
	GAB 30K0.2N	Bench scale, with type approval	30.000	10,00000		
	GAS 3K1DM	Bench scale, with type approval	30.000	0,20000	294 x 225	250,00 €
	GAS 3K1DM	Bench scale, with type approval	1.500; 3.000	1,00000;	230 x 180	180,00 €
	GAS 6K2DM	Bench scale, with type approval	3.000; 6.000	2,00000;	230 x 180	180,00 €
	GAS 15K5DM	Bench scale, with type approval	7.500;	5,00000;	280 x 260	180,00 €
	GAS 30K10DM	Bench scale, with type approval	15.000	10,00000		
	GAS 30K10DM	Bench scale, with type approval	15.000;	10,00000;	280 x 260	180,00 €
	GAS 30K10DM	Bench scale, with type approval	30.000	20,00000		
	GAB-A01	Mains adapter UK	-	-	-	25,00 €
	GAB-A02	Mains adapter US	-	-	-	25,00 €